BORT Generation Rückenorthese

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Generation

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

REF 215 800

Illustrationen Illustrations


1


4


2


5

3

6

Sprachen Languages

DE	Deutsch	Gebrauchsanweisung	04
EN	English	Instructions for use	09
FR	Français	Mode d'emploi	13
ES	Español	Instrucciones de uso	18
IT	Italiano	Istruzioni per l'uso	22
NL	Nederlands	Gebruiksaanwijzing	27
CS	Čeština	Návod k použití	32
ET	Eesti	Kasutusjuhend	36
PL	Polski	Instrukcja użytkowania	40
RO	Română	Instrucțiuni de utilizare	45

 PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for partial relief due to lordosis reduction of the lumbar spine.

Indications

Facet syndrome, lumbosciatica, intervertebral disc protrusions/prolapse, radicular irritation syndrome, degenerative lumbar change, condition following intervertebral disc surgery, stable fractures, spondylosis, spondylolisthesis and spinal canal stenosis.

Contraindications

Unclear soft tissue swellings, impaired sensation and circulation in the specific body region, skin diseases in the specific part of the body.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs

- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

Fitting instructions

Instructions for use for patients:

- 1** When you receive your BORT Generation Back brace, the latch plate is closed. Please open it fully. Also, open both pulley handles and attach them to the side velour surfaces using Velcro.
- 2** Place the back brace around the torso so that the pelvic girdle with its curved protrusion is positioned above the iliac crest. This protrusion is central, on the inside of the pelvic girdle. Moreover, information on the inside of the back brace indicates the direction (**top = towards the head**).
- 3** Close the BORT Generation Back brace on the latch plate applying very little force, using the attached hand straps.
- 4** To close the BORT Generation Back brace fully, grasp the pulley handles and pull them both forwards simultaneously, whereby you only need to apply very little force.

5 Following this, attach the pulley handles to the latch plate or the velour surfaces at the side using the micro Velcro provided. In order to avoid strangulation, the brace should not be closed too tightly. Within the course of this fitting process, the BORT Generation Back brace adapts fully to the anatomy of the body and interlocks.

6 If the lower edge of the brace ends with the hip-bend when sitting, the brace is correctly positioned.

Removal

To remove the brace, loosen the pulley handles off the velour and then open the latch plate.

Info

The back plate is available with a height of 22 cm or 27 cm.

Material composition

polyamide (PA), polyurethane (PU), cotton (CO), elastane (EL), polyester (PES)
The sewn in textile label provides the precise material composition.

Cleaning information

 Hand-washing 30 °C  Do not bleach  Do not dry in a tumble dryer
 Do not iron  Do not clean chemically

Do not use fabric conditioner. Close the Velcro fastening to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly

observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2025

MD Medical device |  Single patient – multiple use

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