	DE Einzelner Patient, mehrfach anwendbar EN Single patient, multiple use FR Patient unique, application polyvalente IT Singolo paziente, uso multiplo PL Jeden pacjent, wielokrotne użycie DA En patient, kan anvendes flere gange
	DE Webseite zur Patienteninformation EN Patient information website FR Site Internet d'informations pour les patients IT Sito web con informazioni per i pazienti PL Strona internetowa z informacjami dla pacjentów DA Hjemmeside til patientoplysninger
	DE Nach Gebrauch wegwerfen EN Discard after use FR Usage unique. Jeter après utilisation IT Smaltire dopo l'uso PL Wyrzucić po użyciu DA Bortskaf efter brug
	DE Verwendbar bis EN Usable until FR Date d'expiration IT Utilizzabile fino a PL Data ważności DA Kan bruges indtil
	DE Hersteller des Produktes EN Manufacturer of the product FR Fabricant du produit IT Produttore dell'articolo PL Producent produktu DA Producent

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INSTRUCTIONS FOR USE EN

INTENDED USE

GENU-TEX® ACL is a soft brace and is used to ensure physiological, two-dimensional guidance of the knee joint and to relieve and protect the joint apparatus in adjustable ranges of motion.

INDICATIONS

The GENU-TEX® ACL can be used for all indications (post-traumatic / post-operative) where it is necessary to secure the physiological guidance of the knee joint, such as:

- Cruciate ligament injuries
- ACL rupture
- Collateral ligament injuries
- General damage to capsule / ligament structures of the knee joint
- Meniscus reconstruction
- Uniaxial or multiaxial instability
- Rheumatoid arthritis
- Knee tep loosening

CONTRA-INDICATIONS

- Allergic, inflammatory or injury-related skin changes (e.g. swellings, redness) of the body areas to be cared for
- Impaired circulation or lymphatic soft tissue swelling
- Neurogenic disorders of the sensory system and skin trophism in the body area to be cared for (sensory disorders with and without skin damage)
- Higher degree of deviation of the leg axis (> approx. 10°)
- Peripheral AVC
- Borderline peripheral blood flow
- Vascular compression risk

SIDE-EFFECTS

If used appropriately and correctly, there are no general side effects, intolerances with longterm effects or allergic reactions known to date.

INTENDED PATIENT TARGET GROUP

The target group is all patients, taking into account indications and contraindications. Instruction and allocation of the correct size is carried out by medical professionals.

PUTTING IT ON AND USE

- To put on the GENU-TEX® ACL knee brace, it is best to sit on a chair so that the knee is in a slight flexion.
- Open all fasteners that have an orange unthreading stop and close them again in the widest setting **1, 2**.
- Take GENU-TEX® ACL by the orange grip aids and step through the orthosis with your foot. This is easier if the back of the orthosis is facing forward. Then pull the brace over the ankle joint to below the hollow of the knee **3**.
- Then turn GENU-TEX® ACL on the lower leg by 180° so that the popliteal knit faces backwards again and the patella pad faces forwards again **4**.

- Then pull the brace over the knee with the help of the grip aid 5 until the silicone ring of the patella pad encloses the kneecap. Check the correct fit of the brace 6 by palpating the kneecap.
- Now close the straps 7. Start with the upper lower leg strap, then follow with the lower thigh strap. Finally, close the knee straps.
- While standing, please check again whether the patella pad continues to enclose the kneecap well. If necessary, correct the position of the brace again and then tighten all the Velcro fasteners so that all four straps fit snugly on the leg with even tension.
- To take off GENU-TEX® ACL, open the front Velcro fasteners 8, close them in the widest setting and slide the brace off the leg 9. It may also help to turn the brace 180° on the lower leg 10.
- Avoid gluing the fasteners to the knitted fabric by always gluing them back onto the belt immediately. This will prevent damage to the knitted fabric.

NOTES FOR THE ORTHOPAEDIC TECHNICIAN

- To adjust the front and rear belt tensions separately, open the Velcro fasteners on the rail pockets.
- The position of the joint splint on the circular belt can now be determined individually.
- The selected position is secured by re-closing the Velcro ears.
- To remove the joint bars, open the proximal closure of joint pocket 1A.
- Loosen the position lock (push button) of the rail 2A.
- The rail can then be removed upwards from the joint pocket 3A, 4A.
- The cover of the joint can be lifted off by rotating it 5A.
- Afterwards, the orange limiting wedges can be removed without tools 6A. The selected wedges can be reinserted in reverse order. By reattaching the joint cover, the wedges are secured in position.
- The joint rail may only be used with inserted wedges. For this purpose, a 0° extension wedge and an F-max flexion wedge are included (delivery condition).
- A lock-out condition is possible in 0°, 10°, 20° and 30°. To do this, first insert the flexion wedge and then the extension wedge.
- If the user has pressure marks in the area of the joint splints, the joint arms can be reshaped in a hot water bath (100° C) from a distance of 3 cm from the joint.
- To adjust the belt tensions separately at the front and rear, the Velcro fasteners at the rail tunnel opened 7A.
- The position of the joint splint on the circular belt can now be individually determined 8A.
- The selected position is secured again by re-closing the Velcro ears 9A.

ACCESSORIES / SPARE PARTS

Optionally available from your specialist store:

- Limitation wedge set ([REF](#) 17720)

IMPORTANT INFORMATION

- Essentially, the indication for and duration of wearing of any orthopaedic device in general, as well as the presence of any of the conditions below in particular, should be discussed with the doctor who is in charge of your treatment.
- Check the functioning of the product with your doctor or specialist retailer.
- The product is designed to treat one patient.
- Do not wear the product directly against open wounds.
- In order for the product to have a long service life and function, it may not be worn in conjunction with oily/greasy or acidic products, ointments or lotions.
- The product can come into contact with splash water.
- Please always close any Velcro fasteners present properly. Damage caused by improper closing of the Velcro fasteners does not constitute a reason for complaint.
- If a significant deterioration in health occurs during the use of the medical device, we ask you to report this serious incident to your specialist dealer, your doctor or us as the manufacturer, as well as to the responsible authority. You can find our contact information in these instructions for use.
- Any improper modifications to the product and/or improper use of the above-mentioned product exclude any product liability on the part of the manufacturer.
- Product may contain metallic objects. We recommend not putting them in contact with electromagnetic interference.
- With products involving direct skin contact with an elastomer (e.g. a silicone pad), blistering may develop on the skin if the product is worn too tightly or if it is moved over the skin with high amounts of friction. In this case, we recommend stopping use of the product. Sweat can develop underneath the silicone, which in rare cases can trigger a reaction to your own sweat with corresponding allergic skin changes. SPORLASTIC is not aware of any allergic reactions being caused by the material. If you have any problems, please talk to your doctor or the supplier of the aid.
- Possible health-related mutual risks or other disadvantages with certain treatments that may arise in association with the use of the product must be discussed with the doctor in charge of treatment.

Please return the packaging and the product to the local recycling collection system.
Please comply with local regulations.

SPORLASTIC GmbH is affiliated with a returns system for packaging and therefore fully complies with the packaging legislation in the Federal Republic of Germany.

CARE

We recommend washing GENU-TEX® ACL gently in lukewarm water (30° C) with mild detergent and air drying. Do not dry the product in the tumble dryer, on the heating or in the microwave. To wash GENU-TEX® ACL, please remove the two joint splints from the brace. Close all Velcro fasteners before washing. This will keep the Velcro fasteners functioning longer.

TRANSPORT AND STORAGE CONDITIONS

Please ensure that the product is stored in a dry place and is protected against moisture and sunlight. Store the product at a normal temperature and humidity.

QUALITY MANAGEMENT SYSTEM

All products from SPORLASTIC GmbH are subject to product inspection as part of our quality management system. If you still have any complaints about our product, we kindly ask that you contact your specialist retailer.

An online version of the instructions for use can be found on our website at www.sporlastic.de.

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