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INSTRUCTIONS FOR USE EN

INTENDED USE

VERTEBRADYN® AKTIV is a back brace for stabilising the lumbar spine with a pressure-distributing back pad.

INDICATIONS

- Conservative / postoperative
- Prolapsed intervertebral disc
- Spinal canal stenosis
- Instability
- Severe lumbar pain
- Muscular imbalances
- Inflammation of the iliosacral joint
- Osteochondrosis
- Facet joint syndrome
- Spondylolisthesis

CONTRA-INDICATIONS

- Allergic, inflammatory or injury-related skin changes (e.g. swelling, redness) of the body areas to be supplied
- Impaired circulation or lymphatic soft tissue swelling
- Neurogenic disorders of the sensory system and skin trophism in the body area to be treated (sensory disorders with or without skin damage)

SIDE-EFFECTS

If used appropriately and correctly, there are no general side effects, intolerances with longterm effects or allergic reactions known to date.

INTENDED PATIENT TARGET GROUP

The target group is all patients, taking into account indications and contraindications. Instruction and allocation of the correct size is carried out by medical professionals.

PUTTING IT ON AND USE

- We recommend wearing a thin, tight-fitting t-shirt under the product.
- Attach the heart-shaped silicone pad in the terry cloth pocket to the middle Velcro velour on side of VERTEBRADYN® AKTIV. The narrow side of the pad should face the orange the instripe on the lower edge of the lumbar brace, with the nub facing the body. The correct positioning of the VERTEBRADYN® AKTIV is achieved when the pad is centred over the lumbar spine 1.
- The front lower end of the VERTEBRADYN® AKTIV should be flush and parallel to the thigh when seated and should not cause any pressure in the groin 2.
- To put on the brace, insert your fingers into the loops on both sides of the VERTEBRADYN® AKTIV closure. Pull the brace apart behind your back with tension and guide it around your body 3.
- Pull your fingers out of the loops and close the Velcro fastener of the VERTEBRADYN® AKTIV completely. When Velcroing, make sure that the closure plates lie positively on top of each other 4.
- Pull the right and left traction element evenly with tension forwards as well as slightly diagonally downwards 5.
- VERTEBRADYN® AKTIV should fit snugly with high compression. In this way, the combination of the stabilising elements in the lumbar region with the pad, the traction elements and the abdominal support can provide significant relief for the spinal column. When fitting with high compression, make sure that it is comfortable to wear.
- Tape the ends of the traction elements under tension to the lower ends of the closure plates 6.

NOTES FOR THE ORTHOPAEDIC TECHNICIAN

- The anatomically pre-shaped spring steel rods can be individually adapted to the patient.

ACCESSORIES / SPARE PARTS

Optionally available from your specialist store:

- BAUCHPLATTE FÜR RÜCKENORTHESEN ([REF] 17428): It can enclose prominent abdominal structures over a large area and thus support the weight of the abdomen.
- VERLÄNGERUNGSTEIL RÜCKEN ([REF] 17400): It is velcroed ventrally to the closure plates and serves as a closure extension for larger abdominal circumferences.
- Silicone pad for VERTEBRADYN® SENSO and VERTEBRADYN® AKTIV ([REF] 27150)
- Bag for silicone pad ([REF] 17150)

⚠ IMPORTANT INFORMATION

- Essentially, the indication for and duration of wearing of any orthopaedic device in general, as well as the presence of any of the conditions below in particular, should be discussed with the doctor who is in charge of your treatment.
- Check the functioning of the product with your doctor or specialist retailer.

- The product is designed to treat one patient.
- Do not wear the product directly against open wounds.
- In order for the product to have a long service life and function, it may not be worn in conjunction with oily/greasy or acidic products, ointments or lotions.
- The product can come into contact with splash water.
- Please always close any Velcro fasteners present properly. Damage caused by improper closing of the Velcro fasteners does not constitute a reason for complaint.
- If a significant deterioration in health occurs during the use of the medical device, we ask you to report this serious incident to your specialist dealer, your doctor or us as the manufacturer, as well as to the responsible authority. You can find our contact information in these instructions for use.
- Any improper modifications to the product and/or improper use of the above-mentioned product exclude any product liability on the part of the manufacturer.
- Product may contain metallic objects. We recommend not putting them in contact with electromagnetic interference.
- With products involving direct skin contact with an elastomer (e.g. a silicone pad), blistering may develop on the skin if the product is worn too tightly or if it is moved over the skin with high amounts of friction. In this case, we recommend stopping use of the product. Sweat can develop underneath the silicone, which in rare cases can trigger a reaction to your own sweat with corresponding allergic skin changes. SPORLASTIC is not aware of any allergic reactions being caused by the material. If you have any problems, please talk to your doctor or the supplier of the aid.
- Possible health-related mutual risks or other disadvantages with certain treatments that may arise in association with the use of the product must be discussed with the doctor in charge of treatment.

Please return the packaging and the product to the local recycling collection system. Please comply with local regulations.

SPORLASTIC GmbH is affiliated with a returns system for packaging and therefore fully complies with the packaging legislation in the Federal Republic of Germany.

CARE

We recommend washing the brace in lukewarm water (30°C) with mild detergent and air drying. Do not dry the product in the tumble dryer, on the heater or in the microwave. To avoid damaging other items of clothing during washing, it is best to wash the brace separately or close the Velcro fasteners beforehand.

Any tension in the knitted fabric that decreases over time is restored by washing. We therefore recommend washing the brace more frequently.

TRANSPORT AND STORAGE CONDITIONS

Please ensure that the product is stored in a dry place and is protected against moisture and sunlight. Store the product at a normal temperature and humidity.

QUALITY MANAGEMENT SYSTEM

All products from SPORLASTIC GmbH are subject to product inspection as part of our quality management system. If you still have any complaints about our product, we kindly ask that you contact your specialist retailer.

MATERIAL COMPOSITION

67 % polyamide, 33 % lycra

An online version of the instructions for use can be found on our website at www.sporlastic.de.

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